



Universal – Flexible - Safe

The DWL® Device Cart was especially designed for the DWL® product line.

Universal and safe!

It is a **universal and safe device cart for examination and treatment rooms, in laboratories, doctors' offices and clinics.**

The **easy-running double steering rolls** easily run over obstacles like cables, hoses and elevator thresholds. The **front rolls can be blocked**, which ensures a safe standing position.

Flexible – maximum space!

The cart offers a **maximum of flexibility and enough space** for device, keyboard, mouse and printer. It also offers a drawer for further accessories.

Due to roll-guided **self-retracting** steel rails the drawers glide easily. They can be taken out for **easy cleaning**.

The optional aluminium pushing/lifting handle facilitates moving and lifting the cart.

The edging of the storage tray ensures a **safe position** of the DWL® devices.

Central ON/OFF switch

6 power supplies can be connected to the rear wall and switched on and off using the **central ON/OFF switch**.

For more detailed information on the DWL® Device Cart, ask the exclusive DWL® distributor for your area or contact us directly on info@dwl.de or by phone on +49 (0) 7731 79769-0.



Advantages of the DWL® Device Cart

- Universal and safe device cart
- Easy-running double steering rolls
- Front rolls can be blocked
- Self-retracting drawer
- Safe position
- 6 power supplies
- Central ON/OFF switch





Article Information

Art.-No.	Description
6933	DWL® Device Cart

Technical Data

Manufacturer	ITD GmbH, Germany (www.itd-cart.com)
Condition at delivery	Completely assembled (delivery of separate parts is not possible)
Load-bearing capacity (on top)	50 kg / 110 lbs
Load-bearing capacity (bottom)	50 kg / 110 lbs
Area (on top)	600 mm x 491 mm (WxD)
Area (bottom)	590 mm x 490 mm (WxD)
Surface	Powder-coated, mouse pad with grey velour surface
Power points	10 A / 250 VAC
Connection cable	Power cable for power points You can select between: <ul style="list-style-type: none"> • Type Cont. Europe <ul style="list-style-type: none"> - power plug: CEE7 - power socket: IEC 320 - nominal current: 10 A - nominal voltage: 250 VAC - according to DIN 40 441 R1/R2 - line H05VV-F 3 G 0.75 mm² (up to a length of 2 m), H05VV-F 3 G 1.0 mm² • Type North America (Hospital Grade) <ul style="list-style-type: none"> - power plug: 3 pins US - socket: IEC 320 - nominal current: 6 A - nominal voltage: 125 VAC - certified according to UL and CSA
Dimensions and weight	840 mm x 705 mm x 625 mm (HxWxD), 38 kg
Dimensions and weight including packaging	1150 mm x 800 mm x 720 mm (HxWxD), 50 kg
Parts delivered	1 top shelf 1 keyboard drawer 1 accessory drawer 1 bottom shelf 1 set of rolls (4 rolls, 2 rolls can be blocked) 1 set of power points 2 power cables 1 connection cable (see above)
Manufactured according to	The DWL carts are manufactured by ITD GmbH according to KD3350.8 and are certified by CE.

Compumedics Germany GmbH
The DWL® Doppler Company
 Josef-Schüttler-Straße 2
 78224 Singen, Deutschland
 t: +49 (0) 7731 79769-0
 f: +49 (0) 7731 79769-99
 e: info@dwl.de
www.dwl.de
www.compumedics.de

Your contact partner for DWL® products

stamp